




# Intubation Aids

Stylets / Introducers / Tube exchanger

# S-Guide

Malleable intubating guide for difficult airway management

## Features




### Soft and Atraumatic Tip

- Preformed
- Coloured:
  - For better visualisation within the glottic structures
  - For easy size recognition



### Oxygenation\* (11 Fr and 15 Fr)

- Three openings at the distal end of the S-Guide to ensure appropriate oxygen flow
- Supplied with two types of O<sub>2</sub> Connectors



### Malleability

- Malleability and rigidity, to preform and manoeuvre the tracheal tube
- Black marking on the sizes 11 Fr and 15 Fr shows ideal position to preform the tracheal tube

Order  
information

**S-Guide** / Two types of O<sub>2</sub> Connectors included (11 Fr and 15 Fr), for single use, sterile

Size	Length	Centimetre Markings	For Tracheal Tube	REF	Box
8Fr	40 cm	5 to 30 cm	≥ 3.0 mm I.D.	<b>33-92-400-1</b>	5
11Fr	65 cm	5 to 40 cm	≥ 4.5 mm I.D.	<b>33-91-650-1</b>	5
15Fr	65 cm	5 to 40 cm	≥ 6.0 mm I.D.	<b>33-90-650-1</b>	5

\*Oxygenation is possible with the sizes 11 Fr and 15 Fr. Oxygenation is not possible with the size 8 Fr.

Additional  
information

Application video



Brochure



Recommendation of use



Annotated bibliography

# / i-Bougie

The versatile Introducer to facilitate difficult intubation

## Features




### Atraumatic Distal Tip

- Rounded and angled design with lateral openings



### Oxygenation

- Allows oxygen flow through the lateral openings at the distal tip
- Supplied with two types of O<sub>2</sub> Connectors



### Shapeability

- Material allows shaping of the i-Bougie
- Centimetre markings from 5 to 40 cm
- Low friction surface for easy railroading of the tracheal tube



#### Order information

**i-Bougie** / Two types of O<sub>2</sub> Connectors included, for single use, sterile

Size	Length	Centimetre Markings	For Tracheal Tube	REF	Box
14Fr	70 cm	5 to 40 cm	≥ 6.0 mm I.D.	33-70-700-1	5

#### Additional information



Application video



Brochure



Annotated bibliography

# X-Changer

Tube exchanger to maintain the airway during tracheal tube exchange




## Features



### Soft and Atraumatic Tip

- Blunt tip is flexible and atraumatic
- Colour coded for easy size recognition



### Oxygenation

- Hollow catheter with lateral openings to deliver oxygen
- Supplied with two types of O<sub>2</sub> Connectors



### Centimetre Markings

- For safe determination of the insertion depth
- Allow proper alignment with the tracheal tube

3 mm



Length 80 cm

VBM

CE 0123

Order  
information**X-Changer** / Two types of O<sub>2</sub> Connectors included, for single use, sterile

Size	Length	Centimetre Markings	For Tracheal Tube	REF	Box
11Fr	80 cm	5 to 40 cm	≥ 4.0 mm I.D.	<b>33-42-800-1</b>	5
14Fr	80 cm	5 to 40 cm	≥ 5.0 mm I.D.	<b>33-41-800-1</b>	5
19Fr	80 cm	5 to 40 cm	≥ 7.0 mm I.D.	<b>33-40-800-1</b>	5

Additional  
information

Application video




Brochure



Annotated bibliography


# / Pocket Introducer

The smallest packaged introducer to facilitate difficult intubation




## Tip

- Preformed and rounded olive tip is atraumatic
- Distal end is shapeable



## Click Connections

- Folded to 20 cm, unfolds to 65 cm with only 3 clicks
- The click connections are permanent and cannot be detached again
- Centimetre markings from 10 to 40 cm



## Smallest Pack Size

- Ideal space solution for emergency bags, helicopters, lab coat pockets, crash carts, etc.
- Ready for use within seconds



## Order information

**Pocket Introducer** / For single use, sterile

Size	Length	Centimetre Markings	For Tracheal Tube	REF	Box
15 Fr	65 cm	10 to 40 cm	≥ 6.0 mm I.D.	<b>33-80-650-1</b>	5

**Pocket Introducer** / In robust plastic cartridge, for single use, sterile

Size	Length	Centimetre Markings	For Tracheal Tube	REF	Box
15 Fr	65 cm	10 to 40 cm	≥ 6.0 mm I.D.	<b>33-81-650-1</b>	1

## Additional information



Application video




Annotated bibliography



Plastic cartridge

# / Stylet

Intubation stylet enables preforming of the tracheal tube



## Features

### Tip

- Soft, atraumatic tip

### Versions

- Without silicone connector
- With moveable silicone connector, which allows fixation of the tracheal tube

Order  
information**Stylet / Reusable**

Size	Length	O.D.	For Tracheal Tube	Without Silicone Connector	Box	With Moveable Silicone Connector (15 mm I.D.)	Box
Small	34 cm	2.0 mm	2.5–3.5 mm I.D.	REF 90-10-111	25	REF 90-10-001	10
Medium	38 cm	3.0 mm	4.0–4.5 mm I.D.	REF 90-10-222	25	REF 90-10-002	10
Large	45 cm	4.0 mm	5.0–6.5 mm I.D.	REF 90-10-333	25	REF 90-10-003	10
Large	45 cm	5.5 mm	7.0–11.0 mm I.D.	REF 90-10-444	25	REF 90-10-004	10

Additional  
information

Annotated bibliography

Additional  
information



Find our full range of products on [www.vbm-medical.de](http://www.vbm-medical.de)

Here you will find our media centre with additional information such as application videos, brochures, flyer etc.

The medical devices in this advertising material are manufactured without the use of natural rubber latex, unless otherwise specified.  
The medical devices in this advertising material do not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.

**VBM** *Medizintechnik GmbH*

Einsteinstrasse 1 / 72172 Sulz a. N. / Germany

Tel.: +49 7454 9596-0 / Fax: +49 7454 9596-33 / e-mail: [sales@vbm-medical.de](mailto:sales@vbm-medical.de) / [www.vbm-medical.de](http://www.vbm-medical.de)

Connect with us:



**CE** 0123

FLY013\_6.0\_EN