

Reusable

Low environmental impact

UDI Code and serial number



Bi-Mask

Sustainable and anatomically shaped anaesthesia face mask


/ Bi-Mask

The Bi-Mask is a face mask made of silicone with integrated plastic shell. This supports a comfortable grip, as the mask is malleable and yet provides stability without collapse. Due to the anatomical shape, the Bi-Mask provides an airtight seal with minimum pressure and the ergonomic design allows an optimised grip via C-grip. The Bi-Mask is available in 5 sizes and processable up to 100 times (autoclavable at 134°C).

Features



Life Cycle Assessment



Life Cycle Assessment (LCA) - September 2025 based on ISO 14040/14044 standards. Centre for Sustainable Healthcare, Amsterdam UMC

Order information

Bi-Mask / Reusable

Size	Patient	Connection	Colour	REF	QTY / Box
2	Child, small	22 mm I.D.	■	35-65-222	1
3	Child, large	22 mm I.D.	■	35-65-223	1
4	Adult, small	22 mm I.D.	■	35-65-224	1
5	Adult, medium	22 mm I.D.	■	35-65-225	1
6	Adult, large	22 mm I.D.	■	35-65-226	1

Fixation Tape / Reusable

Patient	REF	QTY / Box
Child	35-70-155	10
Adult	35-70-255	10

Fixation Tape / For single use

Patient	REF	QTY / Box
Adult	35-30-255	10

Hook Ring / Reusable

REF	QTY / Box
88-60-266	10

This product is manufactured without the use of natural rubber latex, unless otherwise specified.
This product does not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.