

Various openings for a perfect fit of respiration tubing

Tube holder plate swivelling 360°

Made from durable material




## Fixario<sup>®</sup>

Anaesthetic tube holder

# Fixario®

Fixario® is an anaesthetic tube holder for the fixation of respiration tubing, catheters and sampling lines. Its tube holder plate has various openings for a perfect fit of respiration tubing of different size and material, providing strain relief during ventilation and thus relieving pressure on the patient's face.

## Features



## Order information

Fixario® / Anaesthetic tube holder, Length: 49 cm

REF	QTY/Box
73-622-50	1

This product is manufactured without the use of natural rubber latex, unless otherwise specified.  
This product does not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.